

# Filter pipette tips

# Large & Small entrance Available!

- The transparent tip with clear graduations allows easy observation while pipetting.
- Compatible with major pipette brands like Dragon, Brand, Thermo, Eppendorf, and others.
- Vital for liquid manipulation tasks including cell culture and sample prep.











#### ORDER INFORMATION


Size	Packaging	Color	Class	Cat.#	Qty / Case	Size/Case(CM)			Weight/	
3120	Fackaging	Color	Class	Cat.#	Gty / Case	L	W	Н	Case(Kg)	
5mL, Large entrance, Long	24 Tips/Rack	Nature	Common	FTR5000L-S	16 Racks	42.5	28.5	41	6.30	
5mL, Large entrance, Long	100 Tips/Bag	Nature	Common	FTB5000L	10 Bags	44	27	30	3.80	
5mL, Large entrance, Long	24 Tips/Rack	Nature	Ultra Low retention	FTR5000L-L-S	16 Racks	42.5	28.5	41	6.30	
5mL, Large entrance, Long	100 Tips/Bag	Nature	Ollid Low retention	FTB5000L-L	10 Bags	44	27	30	3.80	
5mL, Small entrance	40 Tips/Rack	Nature	Common	FTR5000T-S	16 Racks	42.5	28.5	41	7.20	
5mL, Small entrance	100 Tips/Bag	Nature	Common	FTB5000T	10 Bags	44	27	30	3.80	
5mL, Small entrance	40 Tips/Rack	Nature	Ultra Low retention	FTR5000T-L-S	16 Racks	42.5	28.5	41	7.20	
5mL, Small entrance	100 Tips/Bag	Nature	Ollid Low retention	FTB5000T-L	10 Bags	44	27	30	3.80	
5mL, Large entrance, Short	24 Tips/Rack	Nature	Common	FTR5000E-S	16 Racks	42.5	28.5	41	6.50	
5mL, Large entrance, Short	100 Tips/Bag	Nature	Common	FTB5000E	10 Bags	44	27	30	4.20	
5mL, Large entrance, Short	24 Tips/Rack	Nature	Ultra Low retention	FTR5000E-L-S	16 Racks	42.5	28.5	41	6.50	
5mL, Large entrance, Short	100 Tips/Bag	Nature	Ollid Low retellion	FTB5000E-L	10 Bags	44	27	30	4.20	
10mL, Large entrance	24 Tips/Rack	Nature	Common	FTR10000T-S	16 Racks	42.5	28.5	41	7.10	
10mL, Large entrance	100 Tips/Bag	Nature	Common	FTB10000T	10 Bags	48.5	26	39.5	6.30	
10mL, Large entrance	24 Tips/Rack	Nature	Ultra Low retention	FTR10000T-L-S	16 Racks	42.5	28.5	41	7.10	
10mL, Large entrance	100 Tips/Bag	Nature	Ollid Low retention	FTB10000T-L	10 Bags	48.5	26	39.5	6.30	
10mL, Small entrance	24 Tips/Rack	Nature	Common	FTR10000D-S	16 Racks	42.5	28.5	41	7.00	
10mL, Small entrance	100 Tips/Bag	Nature	Common	FTB10000D	10 Bags	48.5	26	39.5	6.20	
10mL, Small entrance	24 Tips/Rack	Nature	Ultra Low retention	FTR10000D-L-S	16 Racks	42.5	28.5	41	7.00	
10mL, Small entrance	100 Tips/Bag	Nature	Office Low retention	FTB10000D-L	10 Bags	48.5	26	39.5	6.20	
10mL, Small entrance	24 Tips/Rack	Nature	Common	FTR10000E-S	16 Racks	42.5	28.5	41	7.00	
10mL, Small entrance	100 Tips/Bag	Nature	Common	FTB10000E	10 Bags	48.5	26	39.5	6.20	
10mL, Small entrance	24 Tips/Rack	Nature	Ultra Low retention	FTR10000E-L-S	16 Racks	42.5	28.5	41	7.00	
10mL, Small entrance	100 Tips/Bag	Nature	Ollid Low letention	FTB10000E-L	10 Bags	48.5	26	39.5	6.20	


#### COMPATIBILITY

Pipette Brand Cat.#	Eppendorf	Dragon	Sartorius	Brand	Thermo	Watson	Gilson P	Gilson PN	Gilson PL	Gilson PG	Rainin UNV
FTR5000L-S; FTB5000L;	1/	1/	1/								
FTR5000L-L-S; FTB5000L-L	V	V	V								
FTR5000T-S; FTB5000T;				1/	1/						
FTR5000T-L-S; FTB5000T-L				V	V						
FTR5000E-S; FTB5000E;	1/	1/	1/			1/	1/				
FTR5000E-L-S; FTB5000E-L	V	V	V			V	V				
FTR10000T-S; FTB10000T;					1/						
FTR10000T-L-S; FTB10000T-L					V						
FTR10000D-S; FTB10000D;	./	./		./			./				
FTR10000D-L-S; FTB10000D-L	V	V		V			V				
FTR10000E-S; FTB10000E;	1/		1/	1/			1/	1/	1/	1/	1/
FTR10000E-L-S; FTB10000E-L	V		V	V			V	V	V	V	V

# Pipette tips

Large & Small entrance Available!











# ORDER INFORMATION

Size	Packaging	Color	Class	Cat.#	Qty / Case	Size	/Case (0 W	CM) H	Weight/ Case(Kg)
5mL, Large entrance, Long	24 Tips/Rack	Nature		TR5000L-S	16 Racks	42.5	28.5	41	6.20
5mL, Large entrance, Long	100 Tips/Bag	Nature	Common	TB5000L	10 Bags	44	27	30	3.70
5mL, Large entrance, Long	24 Tips/Rack	Nature		TR5000L-L-S	16 Racks	42.5	28.5	41	6.20
5mL, Large entrance, Long	100 Tips/Bag	Nature	Ultra Low retention	TB5000L-L	10 Bags	44	27	30	3.70
5mL, Small entrance	40 Tips/Rack	Nature		TR5000T-S	16 Racks	42.5	28.5	41	7.00
5mL, Small entrance	100 Tips/Bag	Nature	Common	TB5000T	10 Bags	44	27	30	3.70
5mL, Small entrance	40 Tips/Rack	Nature		TR5000T-L-S	16 Racks	42.5	28.5	41	7.00
5mL, Small entrance	100 Tips/Bag	Nature	Ultra Low retention	TB5000T-L	10 Bags	44	27	30	3.70
5mL, Large entrance, Short	24 Tips/Rack	Nature	_	TR5000E-S	16 Racks	42.5	28.5	41	6.40
5mL, Large entrance, Short	100 Tips/Bag	Nature	Common	TB5000E	10 Bags	44	27	30	4.10
5mL, Large entrance, Short	24 Tips/Rack	Nature	***. *	TR5000E-L-S	16 Racks	42.5	28.5	41	6.40
5mL, Large entrance, Short	100 Tips/Bag	Nature	Ultra Low retention	TB5000E-L	10 Bags	44	27	30	4.10
10mL, Large entrance	24 Tips/Rack	Nature	9	TR10000T-S	16 Racks	42.5	28.5	41	7.00
10mL, Large entrance	100 Tips/Bag	Nature	Common	TB10000T	10 Bags	48.5	26	39.5	6.20
10mL, Large entrance	24 Tips/Rack	Nature	TTI T	TR10000T-L-S	16 Racks	42.5	28.5	41	7.00
10mL, Large entrance	100 Tips/Bag	Nature	Ultra Low retention	TB10000T-L	10 Bags	48.5	26	39.5	6.20
10mL, Small entrance	24 Tips/Rack	Nature	G	TR10000D-S	16 Racks	42.5	28.5	41	6.90
10mL, Small entrance	100 Tips/Bag	Nature	Common	TB10000D	10 Bags	48.5	26	39.5	6.10
10mL, Small entrance	24 Tips/Rack	Nature	Ultra Low retention	TR10000D-L-S	16 Racks	42.5	28.5	41	6.90
10mL, Small entrance	100 Tips/Bag	Nature	Ultra Low retention	TB10000D-L	10 Bags	48.5	26	39.5	6.10
10mL, Small entrance	24 Tips/Rack	Nature	Common	TR10000E-S	16 Racks	42.5	28.5	41	6.90
10mL, Small entrance	100 Tips/Bag	Nature	Common	TB10000E	10 Bags	48.5	26	39.5	6.10
10mL, Small entrance	24 Tips/Rack	Nature	Ultra Low retention	TR10000E-L-S	16 Racks	42.5	28.5	41	6.90
10mL, Small entrance	100 Tips/Bag	Nature	Onia Low retention	TB10000E-L	10 Bags	48.5	26	39.5	6.10

# COMPATIBILITY

Pipette Brand											
Cat.#	Eppendorf	Dragon	Sartorius	Brand	Thermo	Watson	Gilson P	Gilson PN	Gilson PL	Gilson PG	Rainin UNV
TR5000L-S; TB5000L;	1/	1/	1/								
TR5000L-L-S; TB5000L-L	V	V	V								
TR5000T-S; TB5000T;				1/	1/						
TR5000T-L-S; TB5000T-L				V	V						
TR5000E-S; TB5000E;	1/	1/	1/			1/	1/				
TR5000E-L-S; TB5000E-L	V	V	V			V	V				
TR10000T-S; TB10000T;					1/						
TR10000T-L-S; TB10000T-L					V						
TR10000D-S; TB10000D;	1/	1/		1/			1/				
TR10000D-L-S; TB10000D-L	V	V		V			V				
TR10000E-S; TB10000E;	1/		1/	1/			1/	1/	1/	1/	1/
TR10000E-L-S; TB10000E-L	V		V	V			V	V	V	V	V