


# L DATA SHEET Tip Combs for Deep Well Plates

## Features



Crafted from high-quality polypropylene, it boasts strong chemical resistance.



## 🖳 Raw Material

These tip combs are produced using pure polypropylene material.

# Application

These products securely position the magnetic rod, enhancing its ability to efficiently adsorb and separate within the sample solution. Polypropylene's inherent low binding affinity towards biomolecules makes it an optimal choice for handling magnetic particles, leading to exceptional recovery of magnetic beads. They are frequently employed alongside deep well plates for nucleic acid extraction experiments.

## We Open To New Ideas!



#### SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

**Add:** Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com

**Web:** www.GenFollower.com **Fax:** +86-575-82032762 **Tel:** +86-575-82936953










# Product Range

CAT	Description	Packaging
DKT0096	96 Tip Combs, sterile, DNase and RNase free	2 PCS/Bag, 5 Bags/Box, 10 Boxes/Case

# Technical Drawing - DKT0096







(Unit: mm)

#### SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

**Add:** Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

**Email:** GenFollower@Genfollower.com

**Web:** www.GenFollower.com **Fax:** +86-575-82032762 **Tel:** +86-575-82936953

