

# DATA SHEET

## Sealing Mat for Deep Well Plates

ISO 9001

### Features



Tailored to accommodate 96-well plates.



The sealing mat fits securely within each well, guarding against the evaporation of volatile samples and cross-contamination.

Sterile and devoid of Dnase, Rnase and human DNA.



### Raw Material

They are made from premium-grade silicone.

### Application

Sealing mats are designed to snugly fit on top of deep well plates, acting as protective shields that cover the plate's wells. This mechanism effectively averts contamination, evaporation, and the potential for cross-contamination. The spectrum of applications for these sealing mats includes safeguarding samples, preventing cross-contamination, ensuring assay uniformity, enabling extended storage, and much more.

**We Open To New Ideas!**



### SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

**Add:** Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

**Email:** GenFollower@Genfollower.com

**Web:** www.GenFollower.com

**Fax:** +86-575-82032762

**Tel:** +86-575-82936953

  
**GenFollower**<sup>®</sup>

# DATA SHEET


## Sealing Mat for Deep Well Plates

ISO 9001

### Product Range

CAT	Description	Packaging
SSM00001	Silicone sealing mat for 96 square well plates	10 PCS/Box, 10 Boxes/Case
SSM00002	Silicone sealing mat for 96 round well plates	10 PCS/Box, 10 Boxes/Case

### Technical Drawing - SSM00001



(Unit: mm)

### SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

**Add:** Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

**Email:** GenFollower@Genfollower.com

**Web:** www.GenFollower.com

**Fax:** +86-575-82032762

**Tel:** +86-575-82936953


**GenFollower**<sup>®</sup>

# DATA SHEET

Sealing Mat for Deep Well Plates

ISO 9001

## Technical Drawing - SSM00002



(Unit: mm)

### SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

**Add:** Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

**Email:** GenFollower@Genfollower.com

**Web:** www.GenFollower.com

**Fax:** +86-575-82032762

**Tel:** +86-575-82936953

 GenFollower<sup>®</sup>