





The wide orifice helps keep delicate samples safe by lowering the force, great for dealing with thick stuff.



Made of polypropylene, it's tough against chemicals.

Compatible with a variety of widely-used pipettes.



The tip's filter effectively prevents cross-contamination.



It's sterilized and clear of Dnase, RNase, DNA, and pyrogens.



Raw Material

The tips are made from medical-grade polypropylene, while the filter in the tip is constructed from polyethylene.

Application

Wide bore tips, also known as large orifice tips, are the go-to choice in labs for tackling thick or dense liquids, like proteins, gels, and genomic DNA. They're especially handy for pipetting high-viscosity samples, working with cell cultures, and handling protein and DNA tasks. These tips are essential for the routine processing of samples, making them a practical choice for a variety of laboratory applications.

We Open To New Ideas!



SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com









Product Range

Volume	Description	Packaging
200uL, Filtered	200uL wide bore filter pipette tips, Nature, sterile, RNase and DNase free	Racked & Bagged
200uL	200uL wide bore pipette tips, Nature, sterile, RNase and DNase free	Racked & Bagged
1000uL, Filtered	1000uL wide bore filter pipette tips, Nature, sterile, RNase and DNase free	Racked & Bagged
1000uL	1000uL wide bore pipette tips, Nature, sterile, RNase and DNase free	Racked & Bagged

Technical Drawing - 200μL, Filtered





SHAOXING GENFOLLOWER BIOTECH CO.,LTD.


Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China


Email: GenFollower@Genfollower.com





Β Technical Drawing - 200μL





(Unit: mm)

SHAOXING GENFOLLOWER BIOTECH CO.,LTD.


Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China


Email: GenFollower@Genfollower.com





Technical Drawing - 1000µL, Filtered





(Unit: mm)

SHAOXING GENFOLLOWER BIOTECH CO.,LTD.


Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China


Email: GenFollower@Genfollower.com





Technical Drawing - 1000μL





(Unit: mm)

SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com

