

Features



The blade, capable of a 360° rotation, can effectively access the growth area, maximizing cell collection.



The soft blade design minimizes potential cell damage.

Four different lengths are offered to accommodate various needs.



Individually packaged units help prevent cross-contamination.



Our product is sterilized and free from DNase, RNase, DNA, heavy metals, and pyrogens.



Raw Material

The handle is crafted from polypropylene, while the blade is produced using thermoplastic elastomers (TPE).

Application

This product is indispensable for the efficient extraction of tissue cultures from flasks, plates, or dishes, resulting in optimal sample recovery. The cell scraper guarantees meticulous collection of intact cell monolayers as it smoothly glides across the growth surface, thereby maintaining the integrity of collected samples throughout the entire harvesting process.

We Open To New Ideas!



SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com









Product Range

CAT.#	Description	Packaging
GD00018	Cell Scrapers, 180mm, Individual packaging, sterile, RNase and DNase free	100 PCS/Bag, 4 Bags/Case
GD00025	Cell Scrapers, 250mm, Individual packaging, sterile, RNase and DNase free	100 PCS/Bag, 2 Bags/Case
GD00235	Cell Scrapers, 235mm, Individual packaging, sterile, RNase and DNase free	100 PCS/Bag, 2 Bags/Case
GD00030	Cell Scrapers, 300mm, Individual packaging, sterile, RNase and DNase free	100 PCS/Bag, 2 Bags/Case

Color







Handle: Yellow/Blue/Green


SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China Email: GenFollower@Genfollower.com





Technical Drawing - GD00018




SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China **Email:** GenFollower@Genfollower.com





Technical Drawing - GD00025



SHAOXING GENFOLLOWER BIOTECH CO.,LTD.


Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com





Technical Drawing - GD00235



SHAOXING GENFOLLOWER BIOTECH CO.,LTD.


Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com





Technical Drawing - GD00030



SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com



L DATA SHEET Cell Lifters







The incorporation of wider blades significantly enhances efficiency.



Designed with straight blades, these lifters effectively minimize the potential for cell damage. Each lifter is individually packaged to ensure no cross-contamination.



The double-blade design caters to distinct needs.



Sterilized and free from pyrogens, DNase, RNase, and DNA.



Raw Material

The lifter is crafted from polypropylene material.

Application

Cell lifters are purposefully crafted tools designed to delicately detach adherent cells from the surfaces of dishes and well plates. This enables the smooth transfer of cells to fresh culture vessels, effectively averting overcrowding and ensuring optimal cell growth conditions. These tools find their primary utility in a diverse range of experiments, encompassing tasks like subculturing, determination, and comprehensive analysis.

We Open To New Ideas!




SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com








Product Range

CAT.#	Description	Packaging
CD00018	Cell Lifters, Single blade, 180mm, Individual packaging, sterile, RNase and DNase free	100 PCS/Bag, 4 Bags/Case
CD00215	Cell Lifters, Double blade, 215mm, Individual packaging, sterile, RNase and DNase free	100 PCS/Bag, 2 Bags/Case

Technical Drawing - CD00018





(Unit: mm)

SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China


Email: GenFollower@Genfollower.com





L DATA SHEET Cell Lifters

Technical Drawing - CD00215



(Unit: mm)

SHAOXING GENFOLLOWER BIOTECH CO.,LTD.

Add: Eastern industrial park, Baiguan street, Shangyu district, Shaoxing city, Zhejiang province, China

Email: GenFollower@Genfollower.com

