


## Centrifuge tubes



- We currently supply 5 types of screw cap centrifuge tubes to meet various needs.
- 1. Achieve a maximum centrifugal force of up to 12,000xg.
- 2. Clear graduations ensure accuraty of samples.
- 3. Highly visible marking areas make it easy to distinguish and classify samples.
- 4. Keeping statement of Dnase, Rnase and pyrogenic free.



## ORDER INFORMATION

Size	Packaging	Туре С	Color of cap	Cat.#	Quantity per Case	Size/Case(CM)			Weight/
						L	W	Н	Case(Kg)
15mL	25 PCS/Bag	Conical	Yellow	B-CTB15YE	20 Bags (25 PCS/Bag, 20 Bags/Case)	41	31	28	4.30
15mL	50 PCS/Foam Rack	Conical	Yellow	B-CTF15YE	10 Racks (50 PCS/Rack, 10 Racks/Case)	44	38	29	4.95
25mL	25 PCS/Bag	Conical	Blue	CTB25BL	8 Bags (25 PCS/Bag8 Bags/Case)	42	35	20.5	2.50
25mL	25 PCS/Foam Rack	Conical	Blue	CTF25BL	8 Racks (25 PCS/Rack\$ Racks/Case)	42	35	20.5	2.80
50mL	25 PCS/Bag	Conical	Yellow	B-CTB50YE	20 Bags (25 PCS/Bag, 20 Bags/Case)	55	32	46	7.75
50mL	25 PCS/Foam Rack	Conical	Yellow	B-CTF50YE	20 Racks (25 PCS/Rack, 20 Racks/Case)	42	34.5	63.5	8.50
50mL	25 PCS/Bag	Self-standing, Conica	l Yellow	B-2SCTB50YE	20 Bags (25 PCS/Bag, 20 Bags/Case)	52	48	29	8.45
50mL	25 PCS/Bag	Self-standing, Conica	l Yellow	B-SCTB50YE	20 Bags (25 PCS/Bag, 20 Bags/Case)	52	48	29	8.45